DK50 DE

INSTALLATION, OPERATION AND MAINTENANCE MANUAL
ИНСТРУКЦИЯ ПО УСТАНОВКЕ, ОБСЛУЖИВАНИЮ И УХОДУ
DIE ANLEITUNG ZUR INSTALLATION, BEDIENUNG UND WARTUNG
NÁVOD NA INŠTALÁCIU, OBSLUHU A ÚDRŽBU

Ⓒ GB

Ⓒ RU

Ⓒ D

Ⓒ SK

CE 0434

**MEDICAL COMPRESSOR
МЕДИЦИНСКИЙ КОМПРЕССОР
MEDIZINALKOMPRESSOR
MEDICINÁLNY KOMPRESOR**

DK50 DE

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07/2012

CONTENTS

1. GENERAL INFORMATION.....	2
PURPOSE.....	2
OPERATOR'S RESPONSIBILITY FOR PATIENT SAFETY	2
MARKINGS	2
WARNINGS.....	2
General safety warnings	3
Electrical system safety warnings	3
WARNING NOTICES AND SYMBOLS	4
USE.....	5
STORAGE AND TRANSPORT	5
2. EQUIPMENT DESCRIPTION.....	6
3. TECHNICAL DATA.....	8
4. OPERATION	9
INSTALLATION AND FIRST OPERATION	9
Removal of transport stabilizers	9
Wheel installation	10
Compressed air connection	10
Electrical connection	10
First operation	11
OPERATION	11
Switching the compressor on	12
Running the compressor	12
Cleaning and replacing filters	12
Cleaning the compressor	12
5. MAINTENANCE.....	13
REPAIRS AND SERVICE	13
Cover removal.....	13
MAINTENANCE SCHEDULE.....	13
Safety valve check	14
Checking tightness of joints and inspecting the equipment.....	15
Replacing filter elements.....	15
Setting output air pressure	16
Cleaning the pressure regulator	16
Replacing the suction filter (18).....	16
Stabilizing the compressor before shipping.....	17
SHUTDOWN	17
EQUIPMENT DISPOSAL	17
6. TROUBLESHOOTING	18
7. SPARE PARTS	19
8. ELECTRIC AND PNEUMATIC DIAGRAMS	20
WIRING DIAGRAM	20
PNEUMATIC DIAGRAM	21
9. PARTS LIST	88

GENERAL INFORMATION

1. GENERAL INFORMATION

PURPOSE

The EKOM DK50 DE is a medical air compressor that supplies clean, oil-free compressed air for use with medical ventilators.

OPERATOR'S RESPONSIBILITY FOR PATIENT SAFETY

The Installation, Operation and Maintenance Manual is an integral part of the equipment and must be kept with the compressor. Careful review of this manual will provide information necessary for correct operation of the equipment.

Rx only

US Federal law restricts the sale of this device by or on the order of a physician.

MARKINGS

Products marked with the CE mark of compliance meet the safety requirements of the European Union (93/42/EEC).

WARNINGS

- The safety of operating personnel and trouble-free operation of the equipment are ensured only if original parts are used. Only accessories and spare parts mentioned in the technical documentation or expressly approved by the manufacturer can be used.
- If any other accessories or consumable materials are used, the manufacturer cannot be held responsible for the safe operation and functionality of the equipment.
- The warranty does not cover damages resulting from the use of accessories or consumable materials other than those recommended by the manufacturer.
- The manufacturer assumes responsibility for the safety, reliability and function of the equipment only if:
 - Installation, calibration, amendments, extensions and repairs are performed by the manufacturer, one of its representatives or a service provider authorized by the manufacturer
 - The equipment is used in accordance with the Installation, Operation and Maintenance Manual
- The Installation, Operation and Maintenance Manual accurately describes the design of the compressor and its compliance with safety and technical standards. The manufacturer reserves all rights to its wiring diagrams, procedures and names.

General safety warnings

The equipment is designed to operate safely when used correctly. Please note the following safety measures to avoid injury or damage.

- Equipment operation must comply with all local codes and regulations.
- Original packaging should be kept for the possible return of the unit. Only original packaging ensures optimal protection of the equipment during transport. If it is necessary to return the equipment during the warranty period, the manufacturer is not liable for damages caused by incorrect packaging.
- The user must immediately notify the supplier if any problem occurs during the use of the equipment.
- This product is not intended for use in areas where there is a risk of an explosion. Do not operate the compressor in the presence of flammable anesthetics.
- Never feed oxygen or nitrous oxide into the compressor. Compressor components are not approved for oxygen or nitrous oxide use.

Electrical system safety warnings

- The equipment must be connected to ground. In order to assure proper grounding, connect the compressor to a receptacle marked "hospital grade."
- Before the compressor is plugged in, make sure that the voltage and frequency of the mains specified on the equipment are the same as the power mains.
- Before operating, check for possible damage to the equipment and any connections. Damaged pneumatic and electrical lines must be replaced immediately.
- If a technical failure occurs, immediately disconnect the equipment from the mains (pull out the main power plug).
- During repairs and maintenance, ensure that:
 - The main power plug is removed from the power socket
 - Compressed air lines are disconnected
 - All pressure has been released from the air tank
- Only a qualified technician can install this equipment.

GENERAL INFORMATION

WARNING NOTICES AND SYMBOLS

The following symbols are used for important information in the Installation, Operation and Maintenance Manual and on packaging and the product:

Attention, see instructions for use

Caution, risk of electric shock

Consult instructions for use

CE mark of compliance

Caution, hot surface

Handling mark on packaging – Fragile, handle with care

Handling mark on packaging – This way up (upright position)

Handling mark on packaging – Protect against moisture

Mark on packaging – Recyclable material

Handling mark on packaging – Temperature during storage and transport

Handling mark on packaging – Limited stacking


Alternating current

Ground connection

Equipotentiality

Fuse

Condensate drain

USE

- The equipment can be installed and operated only in a dry, ventilated and dust-free area. Climatic conditions for operation - see Technical data.
- The compressor must stand on a flat and stable base.
- The compressor must not be exposed to rain. The equipment must not be used in humid or wet environments. Never use the compressor in the presence of flammable liquids or gases.
- Before connecting the compressor to respiration equipment, make sure that it meets the requirements of the respiration equipment. Refer to the Technical data for this purpose.
- Any use other than the compressor's intended use is not considered to be safe. The manufacturer is not responsible for any damages that result if the compressor is used for any other purpose. Risk is exclusively assumed by the operator/user.

STORAGE AND TRANSPORT

The compressor is shipped from the factory in transport packaging with the pump stabilized, protecting it from damage during transport.

For transport, always use the original packaging and secure the compressor in the upright position.

Protect the compressor from humidity, contamination and extreme temperatures during transport and storage. A compressor in its original packaging should be stored in a warm, dry and dust-free area.

Keep the packaging material, if possible. If not, dispose of the packaging material in an environmentally-friendly way. Cardboard can be recycled.

Before moving or transporting the compressor, the pressure in the air tank and hoses must be released and any condensed water must be drained. Secure the motor to prevent movement before shipping.

Prior to transport it is necessary to secure the motor inside the compressor (Chapter 5).

EQUIPMENT DESCRIPTION

2. EQUIPMENT DESCRIPTION

1. Switch, main connector, fuses
2. Hour counter
3. OUT – compressed air output
4. Outlet pressure gauge (accuracy $\pm 5\%$)
5. WALL – input of compressed air from central distribution (Auxiliary equipment)
6. Suction filter
7. Equipotential (ground) pin
8. Condensate tank
9. Power indicator
10. Socket for the electrical cord
11. Compressor
12. Safety valve
13. Air tank
14. Filter (40 μm) and water trap
15. Pressure regulator
16. Filter (5 μm) and water trap
17. Control electronics
18. Suction filter
19. Cooler



The compressor contains an oil-free piston (11) driven by a low-maintenance single-phase electric motor. Compressed air is cooled in the cooler (19) where condensed water is separated into a separate tank (8). Incoming air passes through two filters (6,18) undergoing double filtration as it passes through the system (14,16). Constant pressure at the output is maintained by the pressure regulator (15). The built-in air tank (13) enables peak air consumption of 200 L/min.

The following applies to products with additional equipment – WALL coupling (5):

If the device is equipped with a WALL quick coupling (5), it can be used as a back-up air source. In this configuration, the respiratory equipment is supplied with compressed air from the central air distribution of the medical facility. Air pressure in the central distribution line is sensed by a pressure sensor. If the pressure is sufficient, the compressor stays in STANDBY mode. If the pressure falls, the control unit automatically places the compressor into operation.

When the compressor is used as the main source of air, the control unit determines its operation based on current demand for air. If air consumption is zero, the device switches to STANDBY mode.

EQUIPMENT DESCRIPTION

TECHNICAL DATA

3. TECHNICAL DATA

TYPE	DK50 DE		
VERSION	basic	standard	advanced
Output flow at pressure 3.5 bar (51 psig) L/min	40	50	60
Peak flow	200* L/min (7 Cft/min)		
Voltage / Frequency / Nominal current V/Hz/A	230/50 / 2.8 230/60 / 2.8 120/60 / 5.6 100/50-60 / 5.6	230/50 / 2.8 230/60 / 2.8 120/60 / 5.6 100/50-60 / 5.6	230/50 / 3 230/60 / 3 120/60 / 6 100/50-60 / 6
Air filtration μm	5		
Pressure dew point at 40 L/min (1.4 Cft/min), 3 bar(43.5 psig), 20°C (68°F)	5°C (9°F) below the ambient temperature		
Outlet connection	DISS 1160-A (3/4"-16 UNF)		
Sound level dB(A)	48 49 / 60 Hz	49 50 / 60 Hz	50 51 / 60 Hz
Mode of operation	Continuous - S1		
Separation of condensed water	Automatic		
Output pressure	3.0 bar (43 psig) Adjustable to max. 3.5 bar (51 psig)		
Automatic start up pressure (backup**)	Decrease in pressure in the central distribution line under 2.8 bar (40.6 psig)		
Air tank capacity	2 L (0.61gall UK)		
Pressure range	5 bar (72.5 psig) – 6.5 bar (94 psig)		
Operating pressure of safety valve	7 bar (101.5 psig)		
Adjustment of output air pressure	Pressure regulator		
Dimensions of compressor w x d x h	440 x 350 x 420 mm (17.5 x 14x 16.5 in)		
Dimensions of compressor with wheels w x d x h	470 x 380 x 505 mm (18.5 x 15x 20 in)		
Dimensions of packaging w x d x h	510 x 480 x 470 mm (20 x 19 x 18.5 in)		
Net weight	34 kg (75 lbs)		
Net weight of compressor with wheels	36 kg (79 lbs)		
Gross weight	42 kg (93 lbs)		
Gross weight of compressor with wheels	43 kg (95 lbs)		
Implementation according to EN 60601-1, EN 12021	Type B class I.		
Classification acc. to MDD 93/42 EEC, 2007/47 EC	II b		

* With 0.6 bar drop (8.7 psig)

** Extra equipment

Climatic conditions for storage and transport

Temperature -25°C to +55°C (-13°F to +131°F), 24 hrs +70°C (+158°F)

Relative air humidity 10% to 90% (no condensation)

Relative air humidity for seaworthy packing 10% to 100% (with condensation)

Climatic conditions for operation

Temperature +15°C to +40°C (+59°F to +104°F)

Relative air humidity up to +95%

IPX0 Rating

4. OPERATION

INSTALLATION AND FIRST OPERATION

Do not use the compressor immediately after unpacking as it will not adjust to the ambient temperature.

Only qualified personnel can install the compressor and put it into operation for the first time. The installer shall train the operating personnel in the use and routine maintenance of the equipment. Installation and training of personnel should be acknowledged by the installer's signature on the installation certificate.

Prior to start-up, the four screws for transport stabilization must be removed. If the compressor is switched on without removing the stabilizers, the compressor could be permanently damaged.

The compressor does not contain a backup power supply.

Never block the suction filter on the backside or the vent outlets on the top of the equipment.

If the compressor is equipped with a main source of air, the standby air source must be available.

ANY MODIFICATION OF THIS EQUIPMENT IS FORBIDDEN!


If this equipment is used nearby other instruments, the equipment must be observed in order to verify normal operations in the configuration it will be used.

Instruments may be affected electro-magnetically!

Removal of transport stabilizers

On the bottom of the unit are four M6 stabilizing screws that must be removed before start-up. They are marked with a red warning washer.


Please retain the stabilizing screws for future transport of the compressor.

OPERATION

Wheel installation

Wheels must be installed following the figure below for versions of the compressor with wheels.

Compressed air connection

The given appliance or respiratory apparatus should be connected to the **OUT** (3) quick coupler, as this is the outlet for compressed air.

The following applies for equipment with additional equipment – **WALL** (5) coupling:

Compressed air from the central distribution line is connected via the **WALL** (5) quick coupling as this is the inlet for compressed air, if applicable. Air from the distribution system is automatically connected through the compressor to the **OUT** outlet port. In this configuration, the compressor serves as a backup source of compressed air. If the air pressure from the central distribution system is reduced, the compressor automatically switches on and there is no interruption in the supply of continuously pressurized air at the outlet of the compressor.

Please note that the air supply from the central distribution line entering the compressor must be medical grade air (particulate size, humidity.) The compressor does not modify the air from the central distribution line.

The hose connecting the compressor to the respiration equipment must not pass through a cold environment i.e. placed on the ground. It should be as short as possible with no kinks (this may cause water to condense inside the hose).

Electrical connection

The compressor comes with a plug containing an appropriate protective contact (ground.) Adhere to local electrical regulations. The voltage and frequency of the mains must comply with the specifications on the data label.

The electrical cord must not be stressed or have any tension exerted upon it, and must always be free.

- The socket must be accessible for safety reasons so that the equipment can be safely disconnected from the power supply in case of an emergency.
- The relevant current circuit must be protected.
- Connection of the ground connection (7) to other equipment must adhere to local electrical regulations.
- Fasten the electrical cord through the holder (10).

First operation

- Make sure that the stabilizing screws used during transport were removed.
- Check that the connection to the compressed air supply is correct.
- Check for proper connection to the main power supply.
- Switch on the pressure switch (1) to position "I".

The compressor will be operational once it is put into operation. If there is no air consumption at the output, the device switches to STANDBY mode. When air consumption resumes, the compressor automatically switches on.

During operation, the device drains off the entrapped condensed water from pneumatic circuits via automatic filter separators into a tank.

The following applies to device with additional equipment – WALL (5) coupling:

After the compressor is put into operation, it shall work in one of the following modes depending upon the pressure level in the central distribution and upon air consumption.

- STANDBY – When there is sufficient pressure in the central air distribution line, the main switch is on, the compressor is idle. The device operates as a standby source of air as it checks the pressure in the central air distribution line and, if the pressure falls, the compressor switches on.
- When pressure is low in the central air distribution line, or if the device is not connected to the central distribution line, the compressor is running. If there is no air consumption at the outlet, the device switches to STANDBY mode. When air consumption resumes, the compressor automatically switches on.

During operation, the device drains the trapped condensed water from pneumatic circuits via automatic filter separators into a separate tank.

OPERATION

In case of emergency, switch the equipment off at the switch and pull out the main power plug.

OPERATION


Switching the compressor on

The compressor is switched on at the main power switch (1) by putting it in position "I".

Running the compressor


Compressor operation is automatic. The compressor runs permanently, when air consumption is lower than the set point; once set point is reached the operation is switched to idle run. The pressure gauge displays the value of output pressure. When air consumption stops, the compressor shall switch off after around 30 seconds and switch on after air consumption is restored.

The condensed liquid drains into a tank (8) at the back of the device. When the tank fills up, it must be emptied.

Cleaning and replacing filters

At least once a week take out and clean the suction filter (6) located on the back side. Wash the filter in warm soapy water, rinse thoroughly and allow it to dry. Insert the clean filters so that the intake openings are completely covered by the filters. A replacement suction filter (6) No. 025000018 is included in the basic equipment for the product.

Cleaning the compressor

To clean the compressor, use a detergent that contains no abrasives, chemical solvents or other corrosive agents.

5. MAINTENANCE

REPAIRS AND SERVICE

Warranty and extended warranty repairs are to be completed by the manufacturer or a service provider authorized by the manufacturer.

The manufacturer reserves the right to modify the equipment in any way that will not alter the function or the operation of the equipment.

Only a qualified technician or the Customer Service Department of the manufacturer may perform repairs that go beyond routine maintenance. Use only spare parts and accessories approved by the manufacturer.

Prior to any maintenance or repairs, switch off the compressor and disconnect it from the mains (pull out the main power plug).

Cover removal

- Unscrew the 6 screws from the rear cover
- Disconnect the grounding wire
- Remove the rear cover
- Unscrew the 4 screws from the rear of the main cover
- Disconnect the grounding wire
- Remove the main cover
- Reassemble using the opposite order

MAINTENANCE SCHEDULE


Notice!

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used.

Maintenance	Page	Time interval	To be performed by
Clean suction filter	12	At least once a week	Staff
Safety valve check	14	Once per year	Qualified expert
Replace filters in filter	15	Every 4000 hours	Qualified expert

MAINTENANCE

Test the tightness of joints and inspect the equipment	15	After two years	Qualified expert
Replace piston complete with piston rings, O-rings and bearing	Service documentation	Every 8000 hours	Qualified expert
Replace suction filter	16	Every 8000 hours	Qualified expert
Perform "Repeated Test" according to EN 62353	13	1 x 2 years	Qualified expert


Safety valve check

The safety valve is adjusted to 8 bar by the manufacturer, then tested and marked. It must not be readjusted.

- Turn the screw of the safety valve (12) several rotations to left until the safety valve puffs.
- Only allow the safety valve freely puff **for a short period of time**.
- Turn the screw to the right until it is tight, closing the valve.

Never use the safety valve to release pressure in the air tank. This could damage the safety valve.

Checking tightness of joints and inspecting the equipment

Test for leakage:

- Disconnect the outlet hose from the quick coupling OUT (3); compressed air consumption will stop.
- Pressurize the air tank by releasing some air via a connected device.
- With the compressor off, wait at least 10 minutes before checking the pressure level.
- If the pressure has dropped, use soapy water to find the leaky joint. Start with the joints outside the case, and then test the ones inside the case after removing the covers.
- Tighten or re-seal joints as necessary.

Inspect the equipment:

- Check the status of the compressor motor for balanced operation or noise.
- Check the condition of the hangers above the pump.
- Check fan functionality.
- Check the condition of filters (14,16). With regular operation, condensate drains automatically. Verify this function by comparing the level of condensate in the tank to the scale line indicating the maximum level. If the level of condensate is above the scale line, replace the faulty parts.
- Examine the pump for:
 - Defects in the crankcase
 - Free movement of the crank shaft
- If necessary, replace any faulty parts.

Replacing filter elements

Before beginning, evacuate all compressed air from the tank, reducing the tank pressure to zero and disconnect the equipment from the power mains.

MAINTENANCE

Unscrew the cover and pull it out.
Unscrew the filter holder.
Replace the filter and screw the filter holder into place.
Screw the cover back in place.

Spare filtration pads:

40 µm - white (14)

Item 025200204

5 µm- white-red (16)

Item 025200205


Setting output air pressure

To unlock the control button of the regulator (15), raise it slightly and rotate it to set the demand outlet pressure. After setting the pressure, lock the control button of the regulator by pushing it down.

Cleaning the pressure regulator

- Set the pressure to the minimum position
- Remove the plug (1)
- Rotate and remove the latch (2)
- Remove the line from the valve (3) and the valve (4), then clean the valve seat and the o-ring
- Re-install the components in the following order: (4), (3), (2) and (1)

Replacing the suction filter (18)

- Unlock the cover and pull it out.
- Replace the filter
- Lock the cover


Spare part:

Suction filter 05W POLYESTER, Item 025200194

Stabilizing the compressor before shipping

Before shipping, the compressor must be stabilized to prevent movement. Prior to fixation, it is necessary to disassemble the plastic cover. The air pump must be fixed to the base using the four original M6 screws (see picture). Assemble the plastic cover.

SHUTDOWN

If the compressor is not going to be used, disconnect it from the main power supply and release the pressure in the air tank.

EQUIPMENT DISPOSAL

- Disconnect the equipment from the main power supply.
- Release the pressure in the air tank.
- Dispose of the equipment according to local regulations.
- Parts used in this product have no negative impact on the environment when disposed of properly.

TROUBLESHOOTING

6. TROUBLESHOOTING

Prior to servicing the equipment, reduce the pressure in the air tank to zero and disconnect the equipment from the main power supply.

Only trained service personnel can perform the activities listed in the troubleshooting guide.

PROBLEM	POTENTIAL CAUSE	SOLVING COMMON PROBLEMS
Compressor does not start	There is pressure in the air tank	In STANDBY mode, the compressor starts when pressure decreases
	No main power voltage	Switched off main breaker in distribution system
		Check supply
		Replace a defective fuse
		Loose wire terminal – tighten
	Power cord defective – replace	
	Interrupted winding of motor, damaged thermal protection	Replace motor
Defective capacitor	Replace capacitor	
Jammed piston or other part	Replace damaged parts	
Low compressor performance	Leakage of air from pneumatic distribution system	Inspect pneumatic distribution system – seal loose connection
	Large consumption of compressed air	Do not exceed max. flow (see Technical data)
	Worn out piston rings	Replace worn out piston rings
	Dirty air suction filter	Replace the dirty air suction filter
	Dirty suction filter	Clean or replace the dirty suction filter
Low output compressor pressure	Leakage of air from pneumatic distribution system	Inspect pneumatic distribution system – seal loose connection
	Large consumption of compressed air	Do not exceed max. flow (see Technical data)
	Worn out piston rings	Replace worn out piston rings
	Dirty air suction filter	Replace the dirty air suction filter
	Dirty suction filter	Clean or replace the dirty suction filter
	Dirty filter elements	Replace the dirty elements
Motor has stopped – thermal protection of motor off	Dirty blow down filter	Clean or replace the dirty blow down filter
	Covered vent holes	Uncover vent holes
	Cooling fans don't work	Replace fans
	Unit is hot, unventilated area	Relocate unit
Compressor is noisy (knocking, metal noises)	Damaged piston bearing or bearing in motor	Replace damaged bearing
	Loose (cracked) belt of the air pump hanger	Replace damaged hanger
Water coming out of outlets	Dirty water trap in filter and filter regulator	Clean or replace the water trap
Compressor malfunction	Failure of control unit	Adjust using service software, replace damaged unit


7. SPARE PARTS

• Suction filter (6)			025000018
• Filtration elements (14)			025200204
• Filtration elements (16)			025200205
• Fuse			
version	230V	T6.3A	038100004
	100V, 120V	T10A	038100005
• Insertion	DISS	1160-A	024000261
• suction filter	05W POLYESTER		025200194

ELECTRIC AND PNEUMATIC DIAGRAMS


8. ELECTRIC AND PNEUMATIC DIAGRAMS

WIRING DIAGRAM

- | | |
|----------|----------------|
| PCB | Circuit board |
| EV1, EV2 | Fan |
| YV1 | Solenoid valve |
| M1 | Electric motor |
| Cb1 | Capacitor |
| HL1 | Glow lamp |
| PH1 | Hour counter |

PNEUMATIC DIAGRAM

1. Suction filter
2. Compressor
3. Solenoid valve
4. Return valve
5. Cooler
6. Filter with water trap
7. Pressure regulator
8. Air outlet
9. Air inlet
10. Air tank
11. Safety valve
12. Shuttle valve
13. Pressure gauge

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